

Plasma Thawing Bath

- ❖ COMMON FEATURES
- ❖ Microprocessor based control system
- ❖ Digital LED display of temperature
- ❖ Digital timer for thawing (Optional)
- ❖ Audio-Visual alarm for temperature deviation
- ❖ Audio alarm for timer end
- ❖ Attractively powder coated mild steel exterior
- ❖ Stainless Steel (SS 304) interior with rounded corners for easy cleaning and disinfection
- ❖ Convenient stainless steel basket for holding plasma bags in upright position
- ❖ Resin Bonded Fiber Glass insulation (40 Kg\m3 Density)
- ❖ Mounted on adjustable rubber legs.
- ❖ Tullu noise less pump for deep thawing and circulation of water
- ❖ 1000 watts immersion heater for fast thawing
- ❖ Axial fan for heat exhaust
- ❖ SS 304 Ball Valve (½") for fast Drain
- ❖ 15 Amps, 230 volt, 1 Phase , 50 Hz , Power cord & Plug



PLASMA THAWING BATH

Suitable for rapid thawing 12 frozen plasma bags upto a temperature of 37°C

Plasma Thawing Bath

Model	Bags	Internal Dimension	External Dimension
SAFE thaw	12	(L) 460mm (W) 320mm	(L) 810mm (W)410mm (H) 350mm



2/628, Rapid Nagar, Kunrathur High Road, Gerugambakkam, Chennai - 602 101
Phone: 044 - 23862980 / 81 / 82 Fax: 044-23862984, E-mail: info@cryoscientific.com
website: www.cryoscientific.com

cryo scientific
systems (pvt) Limited